

DNA extractions

Considerations for fungi

Dr Johanna Rhodes


Learning outcomes

- **DNA extractions for NGS**
 - main steps and considerations for fungi
- **Applications**
- **Challenges and limitations**

Who am I?

Hypermutation

- Fungal hypermutators have an elevated mutation rate due to non-synonymous SNPs in DNA mismatch repair system
- Cryptococcus neoformans*
 - Recurrent patient isolates were whole genome sequenced
 - Phylogenetically dissimilar, and large numbers of SNPs and aneuploidy observed
 - Two nonsense mutations in *MSH2*, one each in *MSH5* and *RAD5* (DNA mismatch repair), resulting in premature stop codons in these three genes
 - Hypermutation leads to adaptation of drug resistance
 - Progress of antifungal treatment, and potentially a resistant phenotype.




Rhodes et al. 2017 G3

Centrum voor Infectieziekten
Radboudumc

Zoom interface: Mute, Start Video, Security, Participants (32), Chat, New Share, Pause Share, Stop Share, Annotate, Remote Control, More. A participant has enabled Closed Captioning. Who can see this transcript?


Genomics workflow



Fungal genomics complexities

Considerations

- >1000 fungal species have been genome sequenced
 - >130k bacterial species
- Genome sizes can range from ~2-180 million nucleotides
 - predicted proteomes 2-35 thousand proteins
- DNA extraction difficulties
 - Lysis buffer may not degrade cell wall containing ergosterol
 - beat beating not good for long read sequencing
- Genomes:
 - haploid, diploid, polyploid, aneuploidy



DNA extractions

Considerations for fungi

- Pure or mixed culture, sample or ancient?
- Short or long read technology?
 - Sonication or bead-beating good to break cell wall -> fragmentation
 - Enzymatic digest -> less fragmentation
- Which life cycle stage?
- Personal recommendations (caveats - liquid culture best, and does depend on species....):
 - Illumina: Lucigen MasterPure DNA extraction plus bead beating
 - ONT/PacBio: in-house protocol for filamentous fungi, Epicentre with modifications for yeasts

DNA extractions

Considerations for fungi

- Pure or mixed culture, sample or ancient?
 - What is your question (research/diagnostic)?
 - Pure - liquid culture
 - Sample - host tissue —> dominate sequence reads
 - Extraction from BAL - High Pure PCR Template Prep kit (Roche)
 - Ancient - CTAB method

DNA extractions

Considerations for fungi

- Pure or mixed culture, sample or ancient?
- Short or long read technology?
 - Sonication or bead-beating good to break fungal ergosterol cell wall -> fragmentation
 - short read sequencing
 - Enzymatic digest -> less fragmentation
 - long read sequencing
 - Personnel recommendations for high-molecular weight DNA:
 - MasterPure with large liquid culture start and no bead beating, keep very cold throughout —> great for yeasts (see ‘Rhodes *et al.* Emerging Microbes & Infection’)
 - MasterPure extraction and QIAGEN Uneasy Blood & Tissue kit for clean up —> works for molds (see ‘Hemmings *et al.* Mycopathologia’)

DNA extractions


Considerations for fungi

- Pure or mixed culture, sample or ancient?
- Short or long read technology?
 - Sonication or bead-beating good to break cell wall -> fragmentation
 - Enzymatic digest -> less fragmentation
- Which life cycle stage?
 - spores, hyphae...

DNA extractions

Considerations for fungi - life cycle


- No single extraction method optimal for all species/life stages
- Spores:
 - Thick cell walls, maybe coated, making them resistant to standard lysis —> bead beating, enzymatic treatments
 - Small spore size means less biomass = lower DNA yield
 - Also consider purification to remove contaminants from extracellular material



DNA extractions

Considerations for fungi - life cycle


- No single extraction method optimal for all species/life stages
- Hyphae:
 - Filamentous, less thick cell walls but still a problem —> bead beating
- DNA is usually of high yield and quality/easier to purify



DNA extractions

Considerations for fungi - life cycle


- No single extraction method optimal for all species/life stages
- Mycelium:
 - Dense network of hyphae containing extracellular proteins + polysaccharides which can inhibit extraction —> CTAB or phenol-chloroform
- High biomass = good yield



DNA extractions

Considerations for fungi - life cycle

- No single extraction method optimal for all species/life stages
- Fruiting body
 - Complex, multicellular structure that produces spores.
 - Dense tissue needs extensive mechanical disruption
 - Needs extensive purification e.g. column purification
- Large biomass = good yield



DNA extractions

Considerations for fungi

- Other sequencing approaches:
 - amplicon-based on the MinION (or other ONT platform)
 - QIAGEN PowerSoil Pro kit
 - bead beating using FastPrep-24 5G (Thermo Fisher)

DNA extractions

Other considerations

- Cell wall composition varies across species, life stages and environmental conditions
- Secondary metabolites and pigments
 - can inhibit extractions, so important to perform purification step(s) to improve DNA quality
- DNA yield varies across species
 - Adjust initial sample based on expected/desired yield
- Biosafety requirements
- Storage
- Optimise buffers to improve DNA stability